

2025 VERSION

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments.

www.surgiwelt.com

Implantology
Instruments


SurgiWELT
Medical Instruments Manufacturer



+49 1590 6127425



info@surgiwelt.com



COMPANY PROFILE

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments. With a focus on precision, quality, and innovation.

We provide essential tools for healthcare professionals worldwide. Our products are designed to meet the rigorous demands of modern surgery, ensuring optimal performance and patient safety.

Key areas of expertise:

Manufacturing a wide range of surgical instruments

Commitment to quality and precision

Innovation in surgical instrument design



INDEX

Hand instruments usually double ended, used to separate the periosteum from the underling bone.



1
LaBanca
Periosteal

IMPLAB1
1920

2
LaBanca
Modified

IMPLABM2
1921

Buser

IMPBUS
1922

P9/PA


IMP9/PA
1923

25

IMP25
1924



Hand instruments used to remove hard tissues such as enamel or bone by chipping, cleaving or paring. The beveled blade may be in line with the handle or at an angle to it.



1/2
Buser Modified
Periodontal
Chisel
3mm/4mm

IPCBM1/2
1970

2/3
Buser Modified
Periodontal
Chisel
4mm/5mm

IPCBM2/3
1971

1/3
Buser Modified
Periodontal
Chisel
3mm/5mm

IPCBM1/3
1972

1mm/4mm
Jovanovic
13KTG
Periodontal
Chisel

IPCJOV1/4
1973

2mm/5mm
Jovanovic
13KTG
Periodontal
Chisel

IPCJOV2/5
1974

3mm/6mm
Jovanovic
13KTG
Periodontal
Chisel

IPCJOV3/6
1975

P24G/13KL
Paiacci Modified

IPC24G/13KL
1976

Bone Scopes: Hand instruments similar to a bone chisel but having a hollowed-out blade similar to coupland chisel.

Bone Pluggers: Used to plug and condense grafted bone.

Bone Scopes



Jovaanovic
None Scoop
5mm

IBSJOV5mm
2020

Bone Graft
Scoop

IBSBG
2021

LaBanca
Bone Curette
7mm

ICLBC
2022



LaBanca
3x1.5mm

IBPLAB
2023

4
Jovanovic

IBPJ4
2024

Bone
Packer DE
5 - 7mm

IBPP5-7
2025

Bone Packer
IBP3/5mm
2026

IBPS3-5
2026

Sinus Lift Instruments



Sinus Cobra 3.5mm
ISC3.5
2052



Sinus Cobra In/Out
ISC1/0
2053



Mushroom Sinus Lift
(2-4-5-6-8-10)
IMSL2.8/3.8
2054



Mushroom Sinus Lift
(2-4-5-6-8-10)
IMSL4.8/5.8
2055



Bone Graft (1.0/1.7mm)
(8-11-13-15-18)
IBG1.0/1.7
2081



Bone Graft (1.7/2.0mm)
(8-11-13-15-18)
IBG1.7/2.0
2082



Bone Graft (2.4/2.7mm)
(8-11-13-15-18)
IBG2.4/2.7
2083



Bone Graft (3.3/4.0mm)
(8-11-13-15-18)
IBG3.3/4.0
2084



Bone Graft (3.3/4.0mm)
(8-11-13-15-18), Serrated
IBG3.3/4.0
2085



TPI1-Tunneling Procedure Instrument
ITPI1
2090



TPI2-Tunneling Procedure Instrument
ITPI2
2091



TPI3-Tunneling Procedure Instrument
ITPI3
2092



TPI4-Tunneling Procedure Instrument
ITPI4
2093



TPI5-Tunneling Procedure Instrument
ITPI5
2094



TPI6-Tunneling Procedure Instrument
ITPI6
2095



TPI7-Tunneling Procedure Instrument
ITPI7
2096



Tunneling Knife 1 (Short Handle)
ITK1
2097



Tunneling Knife 1
ITK1
2098



Tunneling Knife 2 (Short Handle)
ITK2
2099



Tunneling Knife 2
ITK2
2100



Giles De Quincey Knife 1
IGDQK1
2101



Giles De Quincey Knife 2
IGDQK2
2102



Steigmann Soft Tissue Splitter
ISSTS
2103



Steigmann Soft Tissue Splitter, Modified
ISSTSM
2104

Osteotomes Chisels: Surgical instruments used to cut or divide bone.

I-8



Straight Osteotome
4mm
IOS4
2110



Curved Osteotome
4mm
IOC4
2111



Straight Osteotome
7.5mm
IOS7.5
2112



Curved Osteotome
7.5mm
IOC7.5
2113



Jovanovic Osteotome
6mm
IOJ6
2114



Surgical/Chisel
3mm S/CUT
IBPPASP
2115



Surgical/Chisel
3mm D/CUT
IBPPASP
2116



OSTM 6mm STR
7-10-13-15-18mm
IBPPASP
2117



OSTM 6mm STR
7-10-13-15-18mm
IBPPASP
2118



OSTM 9mm STR
7-10-13-15-18mm
IBPPASP
2119



OSTM 9mm STR
7-10-13-15-18mm
IBPPASP
2120

Martensite stainless steel tissue punches.

I-9



Tissue Punch CVD
2mm
IBPPASP
2133



Tissue Punch CVD
3mm
IBPPASP
2134



Tissue Punch CVD
4mm
IBPPASP
2135



Tissue Punch CVD
5mm
IBPPASP
2136



Tissue Punch CVD
6mm
IBPPASP
2137



Tissue Punch STR
4mm
IBPPASP
2130



Tissue Punch STR
5mm
IBPPASP
2131



Tissue Punch STR
6mm
IBPPASP
2132

Martensite stainless steel tissue punches.

I-10




Designed to separate the sinus lining from the bone without fraying or tearing.

I-11




Designed to separate the sinus lining from the bone without fraying or tearing.

I-12



Designed to separate the sinus lining from the bone without fraying or tearing.

I-13









CASTROVIEJO
STR 17CM
CVC17CM
2382



Castroviejo
Calipers
Straight
0-40mm
18cm/7"
ICCS1B
2383



Castroviejo
Calipers
Angled
0-40mm
18cm/7"
ICCA18
2385



- 2386-1mm
- 2386-2mm
- 2386-2.5mm
- 2386-3mm
- 2386-4mm
- 2386-5mm
- 2386-6mm
- 2386-7mm (Shown)

NOTE: Also available in 2.5mm, 3mm, 4mm and 6mm

Bone Grafting
Syringe
7mm
IBGS7
2386



Weiss Modified
Castroviejo
Calipers
0-40mm
17cm/6"
ICWCM17
2384



Castroviejo
STR 2.5cm
ICVC-ST8.5cm
2388



Castroviejo
CVD 17cm
ICVC-CV17cm
2389



Cortical
Bone Clamps
IBCC
2387



Bone Grafting 1 Ring Style Syringe
Straight, 3.75mm
2390



Bone Grafting 1 Ring Style Syringe
Straight, 5.0mm
2391



Bone Grafting 1 Ring Style Syringe
Curved, 3.75mm
2392



Bone Grafting 1 Ring Style Syringe
Curved, 5.0mm
2393



Ridge Mapping, Bone Width 0.25mm
Both Side Reading, 14cm
2394



Caliper Double Sided Reading
Scissors Style
2395



Sliding Implant Caliper 94mm,
1.5mm tip diameter
2396



Ruler
15cm/6"
2397



Bone Collector Chisel, Straight
IBCCS
2398



IBCBS
2399 Bone Collector Blade, Straight

Note: Blade Sold Separate



Bone Collector Blade
IBCB
2401

Bone Collector Blade, Curved
IBCCCV
2400

Micro Surgery Endo Forceps, 75mm

IMSEF75
2402

Micro Surgery Endo Forceps, 95mm

IMSEF95
2403

Micro Surgery Endo Forceps, 110mm

IMSEF110
2404

Bone Morselizer

IBM
2405

Bone Crusher
2406

IBM
2406





+49 1590 6127425



info@surgiwelt.com