

2025 VERSION

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments.

www.surgiwelt.com

Orthodontics
Instruments


SurgiWELT
Medical Instruments Manufacturer



+49 1590 6127425



info@surgiwelt.com



COMPANY PROFILE

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments. With a focus on precision, quality, and innovation.

We provide essential tools for healthcare professionals worldwide. Our products are designed to meet the rigorous demands of modern surgery, ensuring optimal performance and patient safety.

Key areas of expertise:

Manufacturing a wide range of surgical instruments

Commitment to quality and precision

Innovation in surgical instrument design



INDEX



HOLE PUNCH PLIERS
Creates Half-Moon Cutouts in Clear Aligners
2426



TEAR DROP PUNCH PLIERS
Designed to create a notch at the
gingival margin of clear aligners
2427



VERTICAL PUNCH PLIER
Vertical Rectangular Shaped Punch Pliers
2429



HORIZONTAL PUNCH PLIER
Horizontal Rectangular Shaped Punch
2428



**PIN & LIGATURE CUTTER
LONG HANDLE**
Cuts Ligature Wires upto 0.015"/0.38mm
2430

15°
Angled



**PIN & LIGATURE CUTTER
NARROW TIP 15, ANGLED**
Cuts Ligature Wires upto 0.015"/0.38mm
2433



**PIN & LIGATURE CUTTER
MINI STYLE**
Cuts Ligature Wires upto 0.015"/0.38mm
2431



**PIN & LIGATURE CUTTER
STANDARD**
Cuts Ligature Wires upto 0.015"/0.38mm
2434



**PIN & LIGATURE CUTTER
NARROW TIP**
Cuts Ligature Wires upto 0.015"/0.38mm
2432

15°
Angled



**PIN & LIGATURE CUTTER
STANDARD 15 ANGLED**
Cuts Wires upto 0.015"/0.38mm
2435



ORTHO HARD WIRE CUTTER
Cuts Hard Orthodontic
Wires upto 0.028"/0.71mm
Precisely Angled Straight Tip
2436



DISTAL END CUTTER
LONG HANDLE CUT-HOLD
Safely Cuts & Holds
Wires upto 0.025"/0.64mm
2440



ORTHO HARD WIRE CUTTER
15, ANGLED
Cuts Hard Orthodontic
Wires upto 0.028"/0.71mm
2437



FLUSH CUT HOLD DISTAL END CUTTER
STANDARD HANDLE
Featuring a Tungsten Carbide Tip Insert
2441



DISTAL END CUTTER
MINI STYLE CUT-HOLD
Safely Cuts & Holds
Wires upto 0.025"/0.64mm
2438



FLUSH CUT HOLD
DISTAL END CUTTER
LONG HANDLE
Featuring a Tungsten Carbide Tip Insert
2442



DISTAL END CUTTER
STANDARD CUT-HOLD
Safely Cuts & Holds
Wires upto 0.025"/0.64mm
2439



HOW PLIERS STRAIGHT
Precisely Serrated Jaws for a Firm Grip
2450



DEBONDING PLIERS STRAIGHT
Precisely Aligned Straight Tip
2451



FORCE MODULE SEPARATING PLIERS
Smooth Reversible Action Control
2452



COON STYLE LIGATURE PLIERS
Long Slim Tips With Reversible Action Control
2453



ADHESIVE REMOVING PLIERS
Double Sided Tungsten Carbide Tip Insert
Comes With Free Pack of Occlusal Pads
2454



WEINGAR PLIERS LONG HANDLE
Designed For Wires upto 0.030"/0.76mm
2455



HOW PLIERS CURVED
Precisely Serrated Jaws for Firm Grip
2457



POSTERIOR BAND REMOVING PLIERS, LONG TIP
Comes with a free pack of occlusal pads
2458



BRACKET REMOVING PLIERS STRAIGHT
Slim edge for an easy grip on bracket
2459



**BRACKET REMOVING PLIERS
ANGLED**
Slim Edge Jaws for an Easy Grip on Bracket
2460



**POSTERIOR BAND REMOVING
PLIERS, STANDARD TIP**
Comes with a free pack of Accusal pads
2461



**DEBONDING
PLIERS, ANGLED**
60° Angled Tip
2462



**WEINGART PLIERS
90° ANGLED**
Designed for wires
upto 0.030"/0.76mm
2463



**WEINGART PLIERS
SLIM**
Designed for wires
upto 0.030"/0.76mm
2464



**WEINGART UTILITY
PLIERS**
Designed for wires
upto 0.030"/0.76mm
2465



**WEINGART UNIVERSAL
45° ANGLED TIP**
Designed for wires
upto 0.030"/0.76mm
2466



**CROWN & BAND
CONTOURING PLIERS**
Featuring a concave and
convex tip design
2467



**CROWN & BAND
CRIMPING PLIERS**
Slim design tips for
precise crimping
2468



**BALL HOOK
CRIMPING PLIERS**
Designed to crimp hooks, posts &
auxiliary stops on arch wire
2469



**POSTERIOR BAND
REMOVING
PLIERS, SHORT TIP**
Comes with a free pack of
accusal pads
2470



BIRD BREAK PLIERS
Featuring a Cone & Pyramid Shaped
Beak Designed for wires
upto 0.028"/0.71mm
2475



BIRD BREAK PLIERS
Featuring a Cone & Pyramid Shaped
Beak Designed for wires
upto 0.028"/0.71mm
2475



TWEED LOOP FORMING PLIERS
Designed for wire upto 0.020"/0.51mm
2477



**LIGHT WIRE PLIERS WITH
3 GROOVES**
Precisely built for wires
upto 0.020"/0.51mm
2478



**LIGHT WIRE PLIERS WITH
CUTTER & 3 GROOVES**
Precisely built for wires
upto 0.020"/0.51mm
2479



**LIGHT WIRE PLIERS
WITH CUTTER**
Precisely built for wires
upto 0.020"/0.51mm
2480



JARABAK PLIERS
Featuring 3 grooves
for a firm wire grip
2481



JARABAK PLIERS WITH CUTTER
Featuring 3 grooves & cutter at base
2482



LINGUAL ARCH FORMING PLIERS
Designed to create double back and
triple back bends in all wires
upto 0.036"/0.91mm
2483



NANCE LOOP FORMING PLIERS
For step beak to create 3,4,5 and 6mm loop, designed for all wires upto 0.022"/0.56mm
2484



DE LA ROSA PLIERS
Featuring a chamfered concave & convex tip built for all wires upto 0.025"/0.64mm
2485



DE LA ROSA PLIERS WITH 3 GROOVES
Featuring a chamfered concave & convex tip built for all wires upto 0.025"/0.64mm
2486



HOLLOW CHOP CONTOURING PLIERS
Long Concave & Convex tip for intricate wire contouring
2487



OPTICAL PLIERS
Featuring a slim beak with concave end designed for all wire upto 0.025"/0.64mm
2488



V-STOP PLIERS
Create V-Shaped bends on wires upto 0.022"/0.56mm
2489



YOUNG PLIERS
3 Sectional beak with 3 grooves designed for wires upto 0.028"/0.71mm
2490



TWEED ARCH FORMING PLIERS 1.75mm
Designed for wires upto 0.025"/0.66mm
2491



TWEED ARCH FORMING PLIERS 1.0mm
Designed for wires upto 0.025"/0.66mm
2492



**TORQUING PLIERS
FEMALE/OUTER**

Places torques in all wires upto 0.030"/0.76mm
alongside our inner torquing pliers
2493



**TORQUING PLIERS
MALE/INNER**

Places torques in all wires upto 0.030"/0.76mm
alongside our outer torquing pliers
2494



STEP PLIERS 1MM
Designed for all wires
upto 0.028"/0.71mm
2495



STEP PLIERS 0.75MM
Designed for all wires upto 0.028"/0.71mm
2496



STEP PLIERS 0.5MM
Designed for all wires upto 0.028"/0.71mm
2497



**NITI DISTAL
CINCH BACK PLIERS**
Designed for SS & NiTi Wires
upto 0.025"/0.64mm
2498



**NITI DISTAL CINCH BACK
PLIERS, LINGUAL**
Designed for SS & NiTi Wires
upto 0.018"/0.46mm
2499



**3 PRONG PLIERS
STANDARD**
Precisely built for all wires
upto 0.030"/0.76mm
2500



**3 PRONG PLIERS
HEAVY**
Precisely built for all wires
upto 0.030"/0.76mm
2501



**3 PRONG PLIERS
NiTi**
Precisely built for all wires
upto 0.020"/0.5mm
2502



Ligature
Director/Band Pusher
2506



Small
Plugger / Picker
2507



Dull Scalor /
Band Pusher
2508



Ligature
Director/Scaler
2509



Scaler /
Band Pusher
2510



Lingual Ligature
Director U Shape
2511



Dull
Scaler

2512



Dull
Scaler

2513



Large
Plugger/Picker

2514



Picker/
Ligature Director

2515



Scaler /
Picker

2516



Dull
Scaler

2517



Ligature Director/
Band Pusher

2518



Band Pusher
With Short Tip
2519



Band Pusher
2520



Filling Explorer/
Ligature Director
2521



Band Pusher/
Scaler
2522




Band Pusher/
Scaler
2523



Ligature
Director
2524



Ligature
Director
2525



Ligature
Torquing Key
2526

Straight .020"
Applicator
2527

Single End
Ligature Director
2528

Lingual .018"
Ligature Director
2529

Lingual .025"
Ligature Director
2530

Ligature
Director/Pluggers
2531

Ligature
Director/Tucker
2532

Plugger
Ligature Picker
2533

Posterior Bracket Holder Tweezer, Long



2566

Direct Bond Bracket Tweezer With Aligner & Positioner



2567

Buccal Tube Tweezer



2568

Slim Buccal Tube Bonder Tweezer



2569

Direct Bond Bracket Tweezer, Lingual



2570

Direct Bond Tweezer, Aligner



2571



All Purpose Scissors
2572






Crown Scissors
2573



Crown Scissors, Straight
2574



Crown Scissors, Curved
2575





MATHIEU NARROW TIP
Enhanced double spring function
2585



MATHIEU SMAHA FINE TIP
Enhanced double spring function
2586



MATHIEU
NARROW TIP, TC
Enhanced double spring
function
2587



MATHIEU
STANDARD TIP
Enhanced double spring
function
2588



MATHIEU
STANDARD TIP, TC
Enhanced double spring
function
2589



+49 1590 6127425



info@surgiwelt.com