

# 2025 VERSION

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments.

[www.surgiwelt.com](http://www.surgiwelt.com)

**GENERAL SURGICAL**  
*Instruments*



**Surgiwelt**  
*Medical Instruments Manufacturer*



+49 1590 6127425



[info@surgiwelt.com](mailto:info@surgiwelt.com)



## **COMPANY PROFILE**

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments. With a focus on precision, quality, and innovation.

We provide essential tools for healthcare professionals worldwide. Our products are designed to meet the rigorous demands of modern surgery, ensuring optimal performance and patient safety.

Key areas of expertise:

Manufacturing a wide range of surgical instruments

Commitment to quality and precision

Innovation in surgical instrument design



## INDEX

These instruments are used to loosen the tooth from peridental ligament and ease extraction

2.5mm



EZ Luxators  
Small, Straight 2.5mm  
SHLSS2.5mm  
4200

2.5mm



EZ Luxators  
Small, Curved 2.5mm  
SHLSC2.5mm  
4201

2.5mm



EZ Luxators  
Spade-Straight  
SHLSS  
4202

2.5mm



EZ Luxators  
Small, Distal 2.5mm  
SHLSD2.5mm  
4203

4.5mm



EZ Luxators  
Small, Mesial 2.5mm  
SHLSM2.5mm  
4204

4.5mm



EZ Luxators  
Large, Distal 4.5mm  
SHLSD4.5mm  
4205



EZ Luxators  
Large, Mesial 4.5mm  
SHLSM4.5mm  
4206





H5-Hourigan  
Serrated  
SH5  
4215

H6-Hourigan  
Serrated  
SH6  
4216



11 Straight Bernard  
SBERS11  
4207



6 Right Bernard  
SBERR6  
4208



7 Left Bernard  
SBERL7  
4209



Bernard Spade, Offset IN  
SBERO-IN  
4210



Bernard Spade, Offset OUT  
SBERO-OUT  
4211



These instruments are used to loosen the tooth from  
peridental ligament and ease extraction



8  
Crane Pick  
**E8**  
4220



12  
**E12**  
4221



11M  
MacMillan  
Gouge  
**E31F**  
4223

11M  
MacMillan  
Gouge  
**E12M**  
4224

13M  
MacMillan  
Gouge  
**E13M**  
4225



19  
**E19**  
4228



20  
**E20**  
4229



23  
**E23**  
4222



31F  
Friedman  
[Craig]  
Gouge Serrated  
**E31F**  
4226

32F  
Friedman  
[Craig]  
Gouge Serrated  
**E32F**  
4227



38  
Serrated  
**E38**  
4230



39  
Serrated  
**E39**  
4231



44C  
Bradley  
**E44C**  
4232

These instruments are used to loosen the tooth from peridental ligament and ease extraction



46

**E46**

4233



46R

**E46R**

4234



69W  
Woodward Hand  
Pressure Chisel  
**E46R**  
4235



77

**E77**

4236



77R

**E77R**

4237



80

**E80**

4238



81

**E81**

4239



81A

**E81A**

4240



81E

**E81E**

4241



92  
Serreted

**E92**

4242



93  
Serreted

**E93**

4243



98

**E98**

4244



190

**E190**

4245



191

**E191**

4246

These instruments are used to loosen the tooth from  
peridental ligament and ease extraction

Cogswell



A  
Cogswell  
**EA**  
4247

B  
Cogswell  
**EB**  
4248

Stout



11  
Stout  
**E11**  
4249

11A  
Stout  
**E11A**  
4250

Heidbrik



1  
Heidbrik  
Straight  
**E1**  
4251

2  
Heidbrik  
Straight  
**E2**  
4252



40  
Heidbrik  
Straight  
**E40**  
4253

41  
Heidbrik  
Straight  
**E41**  
4254

These instruments are used to loosen the tooth from peridental ligament and ease extraction

### Cryer:

Cryer's right and left elevators with triangular sharp-pointed blades and hollow handles



51  
Curtis

**E51**  
4255

52  
Curtis

**E52**  
4256



44  
Cryer Small

**E44**  
4257



45  
Cryer

**E45**  
4258



25  
Cryer  
Medium

**E25**  
4259



26  
Cryer  
Medium

**E26**  
4260



A30  
Cryer

**EA30**  
4261



A31  
Cryer

**EA31**  
4262



31  
Cryer Large

**EA31**  
4263



32  
Cryer Large

**EA32**  
4264



These instruments are used to loosen the tooth from peridental ligament and ease extraction

**Luxating Elevator:** The luxating elevators have thin, sharp blades that are used to cut periodontal ligaments.

These instruments are used by rocking and cutting verses and elevator that lifts and pries. The importance of luxating elevator is to reduce pressure on adjacent tooth.

### LUXATOR ELEVATORS



Luxating Elevator  
Straight, 2mm  
SLES2mm  
4290



Luxating Elevator  
Curved, 2mm  
SLEC2mm  
4291



Luxating Elevator  
Straight, 2.5mm  
SLES2.5mm  
4292



Luxating Elevator  
Straight, 2.5mm  
SLES2.5mm  
4293



Luxating Elevator  
Straight, 3mm  
SLES3mm  
4294



Luxating Elevator  
Curved, 3mm  
SLEC3mm  
4295



Luxating Elevator  
Straight, 4mm  
SLES4mm  
4296



Luxating Elevator  
Curved, 4mm  
SLEC4mm  
4297



Luxating Elevator  
Straight, 5mm  
SLES5mm  
4298



Luxating Elevator  
Straight, 5mm  
SLES5mm  
4299

### Coupland Gouge



1C  
Coupland  
Gouge  
**E1C**  
4300-1



2C  
Coupland  
Gouge  
**E2C**  
4300-2



3C  
Coupland  
Gouge  
**E3C**  
4300-3



A1C  
Coupland  
Gouge  
**EA1C**  
4301



A2C  
Coupland  
Gouge  
**EA2C**  
4302



A3C  
Coupland  
Gouge  
**EA3C**  
4303



These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

### T-bar Handle

The Hourigan T-bar Handle is twice  
The length of the potts T-bar handle



### Cross-bar Handle



These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

## Heidbrink



## Apical:

Set of three Elevators.. Right, Left and Straight



These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

### West



### Howard

Set of three Elevators.. Right, Left and Straight




These instruments are used to loosen the tooth from the periodontal ligament and ease extraction.

Davis Teasers





Thin sharp blades used to facilitate the removal of teeth with minimal damage to the surrounding alveolar bone by severing the periodontal ligament. May be used wherever an extraction is indicated, but are particularly important when considering the placement of a dental implant.



1 Posterior  
SPT1  
4347

3 Posterior-Right  
SPTR3  
4348

4 Posterior-Left  
SPTL4  
4349

2 Anterior  
SAT2  
4350

5 Anterior-Angled  
SATA5  
4351

6 Anterior-Straight  
SATS6  
4352

Wise  
SPW  
4353



These instruments are used for retracting periosteum after removing gingival tissue.

**S-13**



These instruments are used for retracting periosteum after removing gingival tissue.



3/9A



4409



PR-3

4410



Tin Coating

PR-3

4411



P10

4412



P15

4413



Howard-9

4414



Howard-9

Mini

4415



Allen-9

Mini

4416

These instruments are used for retracting periosteum after removing gingival tissue.





These instruments are used for retracting periosteum after removing gingival tissue.



HF-20

4425



Tin Coating  
HF-20

4426



Molt-5

4427



Selden-23

4428



Freer

4429



OHL-10  
Professional  
Elevator

4430



Woodson-1  
Modify

4431

These instruments are used for curettage and cyst removal and debridement of tooth socket.

**S-17**

## Lucas-Curettes






These instruments are used for curettage and cyst removal and debridement of tooth socket.

**S-18**

## Molt-Curettes


## Moller-Curettes



## BONE CHISELS SINGLE BEVEL



## DOUBLE BEVEL



Bone Files: For final smoothing of bone by pushing out.

### Miller Collburn-Cross Cut



The small paddle shaped end is half the size of the larger paddle shaped end.

1X  
Miller-Collburn  
Cross-Cut

4600



Curved working end combined with a straight working end.

2X  
Miller-Collburn  
Cross-Cut

4601



Curved working end combined with a small paddle shaped working end.

3X  
Miller-Collburn  
Cross-Cut

4602

### Miller-Straight Cut



Same design as the 1X Cross-Cut but with straight cut.

21  
Miller

4603



Same design as the 2X Cross-Cut but with straight cut.

45  
Miller

4604

Bone Files: For final smoothing of bone by pushing out.

### Miller-Straight Cut



52  
Miller

4605



64  
Miller

4606

### Wahl-Straight Cut



1  
Wahl

4607



2  
Wahl

4608

### Howard-Straight Cut



12  
Howard

4609




67  
Howard

4610




Used for trimming tissue or cutting sutures.  
All Scissors Shown 50% Size

**S-24**



1  
Kelly  
Straight/Pointer  
16cm/6 1/4"  
**S1**  
4723



2  
Kelly  
Straight/Pointer  
16cm/6 1/4"  
**S2**  
4724



Kelly  
Curved  
18cm/7"  
**S1L**  
4725



Kelly  
Straight  
18cm/7"  
**S2L**  
4726



5  
Wagner  
Straight  
11.5cm/4 1/2"  
**S5**  
4720





5  
Wagner  
Curved  
11.5cm/4 1/2"  
**S6**  
4721



7  
Wagner  
Angled  
11.5cm/4 1/2"  
**S7**  
4722

Used for trimming tissue or cutting sutures.  
All Scissors Shown 50% Size

**S-25**



Used for trimming tissue or cutting sutures.  
All Scissors Shown 50% Size

**S-26**



31  
Castroviejo  
Curved  
10cm/4"  
**S31**  
4740



17  
Iris  
Straight/Delicate  
11.5cm/4 1/2"  
**S17**  
4741



18  
Iris  
Curved/Delicate  
11.5cm/4 1/2"  
**S18**  
4742



18  
Metzenbaum  
Baby  
12.5cm/5"  
**SMB**  
4743



Joseph  
Curved  
14cm/5 1/2"  
**SJC**  
4735



Joseph  
Straight  
14cm/5 1/2"  
**SJS**  
4736



3  
Mayo  
Curved/Blunt  
14.5cm/5 3/4"  
**S3**  
4737



4  
Mayo  
Curved/Blunt  
14.5cm/5 3/4"  
**S4**  
4738



Strabismus  
Straight  
Tenotomy  
11.5cm/4 1/2"  
**S40**  
4739

Used for trimming tissue or cutting sutures.

All Scissors Shown 50% Size

**S-27**



21  
Straight/  
Pointed  
14.5cm/5.3 " "  
**S21**  
4744



22  
Straight/  
Blunt-Sharp  
14.5cm/5.3 " "  
**S22**  
4745



23  
Curved/  
Pointed  
14.5cm/5.3 " "  
**S23**  
4746



Lister  
Bandage  
14.5cm/5.3 " "  
**SLB**  
4747



13  
Suture  
15cm/6 " "  
**S13**  
4748



13S  
Suture  
15cm/4.3 " "  
**S13S**  
4749



Crown &  
Gold  
Straight  
15cm/4.3 " "  
**SCGS**  
4750



Crown &  
Gold  
Curved  
15cm/4.3 " "  
**SCGC**  
4751



Crown &  
Gold Curved  
Pedodontic  
15cm/3.5 " "  
**SCGCP**  
4752



**Super-Cut Scissors:** The razor-sharp cutting edge cuts effortlessly and softly also through close material. The serrated cutting edge on the second blade prevents a slipping from tissue and sutures. Super-Cut is not practical for wires, cartilages and other hard material. As identification, one handle ring is gold plated.

All Scissors Shown 50% Size

Super Cut Scissors



10  
Double-Curved  
13.5cm/5 1/2"  
**S10SC**  
4800



14  
Double-Curved  
11.5cm/4 1/2"  
**S145C**  
4801



16  
Gold man-Fox  
15.5cm/6 1/8"  
**S165C**  
4802



17  
Iris  
Straight-Delicate  
11.5cm/4 1/2"  
**S17SC**  
4803



18  
Iris  
Curved-Delicate  
11.5cm/4 1/2"  
**S18SC**  
4804



Metzenbaum  
Straight/Blunt  
14.5cm/5 3/4"  
**S5054**  
4805



Metzenbaum  
Curved/Blunt  
14.5cm/5 3/4"  
**S5055**  
4806



Metzenbaum  
Straight/  
Pointed  
14.5cm/5 3/4"  
**S5056**  
4807



Metzenbaum  
Curved/  
Pointed  
14.5cm/5 3/4"  
**S5057**  
4808



Metzenbaum  
Delicate-  
Curved/Blunt  
14.5cm/5 3/4"  
**S5069**  
4809

All Scissors Shown 50% Size

**S-29**



Goldman-Fox  
Straight  
13cm/5"

**S5080**  
4810



Goldman-Fox  
Curved  
13cm/5"

**S5081**  
4811



Wire Cutting  
Notched  
12cm/4"

**S5085**  
4812



Iris (17)  
Straight  
11.5cm/4"

**S5082**  
4813



Iris (17)  
Curved  
11.5cm/4"

**S5083**  
4814



Dean  
16.5cm/6"

**S5009**  
4815



Locklin  
Curved  
Handle  
16.5cm/6"

**S5011**  
4816



Locklin  
Straight  
Handle  
16.5cm/6"

**S5012**  
4817



Kelly  
Curved  
16.5cm/6"

**S5001**  
4818



Kelly  
Straight  
16.5cm/6"

**S5002**  
4819



These pliers are used to grasp and stabilize soft tissue flaps during suturing.

All tissue pliers and forceps are Shown 50% Size

**S-30**

## Semkin-Taylor



31  
Semkin-Taylor  
Straight  
12.5cm/5"  
**SLB**  
4870



32  
Semkin-Taylor  
Curved  
12.5cm/5"  
**SLB**  
4871



33  
Semkin-Taylor  
1x2 Straight  
12.5cm/5"  
**SLB**  
4872



34  
Semkin-Taylor  
1x2 Curved  
12.5cm/5"  
**SLB**  
4873



18  
3X4  
15cm/6"  
**SLB**  
4874



50-L  
Fox Locking 3X3  
15cm/6"  
**SLB**  
4875



Kramer-Nevins  
Tissue Pliers  
1X2 Curved  
**SLB**  
4876

## Suture Pliers



20 Corn  
Suture Pliers  
15cm/6"  
**SLB**  
4877



41  
Adson Plain  
12cm/4 3/8 "  
**SLB**  
4878



42  
Adson 1X2  
12cm/4 3/8 "  
**SLB**  
4879

43  
Adson-Brown  
7X7  
12cm/4 3/8 "  
**SLB**  
4880



Russian  
15cm/6"  
**SLB**  
4881



Russian  
20cm/8"  
**SLB**  
4882



Gerald Plain  
18cm/7"  
Straight/Serrated  
**SLB**  
4883



Gerald 1X2  
15cm/7"  
Straight  
**SLB**  
4884

## Alison



Allison  
Baby 4X5  
12.5cm/5"  
**SLB**  
4885



Allison  
Straight 5X6  
15cm/6"  
**SLB**  
4886

All tissue pliers and forceps are Shown 50% Size

**S-31**

## Adson




Adson  
12cm/4<sup>3</sup>/<sub>8</sub> "  
**TCTPA**  
4887

Adson 1X2  
12cm/4<sup>3</sup>/<sub>8</sub> "  
**TCTPA**  
4888

Adson  
Brown  
12cm/4<sup>3</sup>/<sub>8</sub> "  
**TCTPAB**  
4889

## General



General  
15cm/6"  
**SLB**  
4890

General  
1X2  
18cm/8"  
**SLB**  
4890

## Gerald



Gerald  
18cm/7"  
**SLB**  
4892

Gerald 1X2  
18cm/7"  
**SLB**  
4893

Gushing  
18cm/7"  
**SLB**  
4894

## Cushing



Gushing  
Angled  
18cm/7"  
**SLB**  
4895

Gushing  
Bayone  
19cm/7<sup>1</sup>/<sub>2</sub> "  
**SLB**  
4896

Potts Smith  
18cm/7"  
(Shown)  
**SLB**  
4897

Lister  
Bandage  
14.5cm/5<sup>3</sup>/<sub>8</sub> "  
**SLB**  
4898

Da Bakey  
15cm/6"  
(Shown)  
**SLB**  
4899

De Bakey  
19cm/7<sup>1</sup>/<sub>2</sub> "  
**SLB**  
4900





Martensite stainless steel instruments are used for grasp and guide the needle with white suturing  
All Needle Holders are Shown 50% Size



Used for grasp and guide the needle with white suturing  
All Needle Holders are Shown 50% Size

Halsey



Hegar-  
Baumgartner  
11.5cm/4 1/2 "

**NH5034**  
4961



Halsey  
Smooth  
13cm/5 "

**NH5036**  
4962

Halsey  
13cm/5 "

**NH5037**  
4963



Castroviejo  
Straight  
14cm/5 1/2 "

**NH5020**  
4964



Castroviejo  
Curved  
14cm/5 1/2 "

**NH5021**  
4965



Micro  
Castroviejo  
Straight  
14cm/5 1/2 "

**NH5020M**  
4966



Micro  
Castroviejo  
Curved  
14cm/5 1/2 "

**NH5021M**  
4967



Castroviejo  
Straight  
18cm/7 "

**NH5024**  
4968




Castroviejo  
Straight  
Round Handle  
14cm/5 1/2 "

**NH5024R**  
4969

All Needle Holders are Shown 50% Size

**S-34**



Arruga  
Serrated  
16cm/6" "

**NH5022**  
4970



Webster  
Smooth  
11.5cm/4" "

**NH5030**  
4971



Chile-Wood  
15cm/6" "

**NH5038**  
4972



Mayo-Hegar  
16cm/6" "

**NH5042**  
4976



Olsen-Hegar  
Needle Holder  
17cm/6" "

**NH5068**  
4973



Olsen-Hegar  
Needle Holder  
Scissors  
14cm/5" "

**NH5072**  
4974



De Bakkey  
18cm/7" "  
Delicate

**NH5046**  
4975



Mayo-Hegar  
18cm/7" "

**NH5044**  
4977



Used for grasp and guide the needle with white suturing  
All Needle Holders are Shown 50% Size



Matheu-Kocher  
14cm/5 1/2"

**NH5074**  
4978



Matheu-  
14cm/5 1/2"

**NH5076**  
4979



HMT  
Thin Jaw  
14cm/5 1/2"

**NH5050**  
4980



Fine Swedish  
Pattern  
15cm/6"

**NH5052**  
4981



Micro-Vascular  
Pattern  
15cm/6"

**NH5082**  
4982



Micro-Vascular  
Pattern  
15cm/7"

**NH5084**  
4983



Mini-Ryder  
Pattern  
15cm/6"

**NH5092**  
4984



Mini-Ryder  
Pattern  
18cm/7"

**NH5094**  
4985



A multipurpose instruments used to clamp off blood vessels, remove small root tips and grasp loose objects.  
All Forceps are Shown 50% Size

**S-36**

## Hartman



Hartman  
Straight  
10cm/4"

**HHS**  
5030



Hartman  
Curved  
10cm/4"

**HMC**  
5031



5  
1x2 Teeth  
Straight  
12.5cm/5"

**HHS**  
5034

## Adson



Adson  
Straight  
19cm/7"

**HAS**  
5032



Adson  
Curved  
19cm/7"

**HAC**  
5033



Rochester-Pean  
Curved  
16cm/6"

**HRPC**  
5035

A multipurpose instruments used to clamp off blood vessels, remove small root tips and grasp loose objects.  
All Forceps are Shown 50% Size

**S-37**

## Halsted-Mosquito



3  
Halsted-Mosquito  
Straight  
12cm/4"  
**HHS**  
5036



3  
Halsted-Mosquito  
Curved  
12cm/4"  
**HHC**  
5037



Kocher  
1x2 Teeth Curved  
14cm/5"  
**HKC**  
5038

## Kelly



1 Kelly  
14cm/5"  
**HHS**  
5039



2 Kelly  
Curved  
14cm/5"  
**HHS**  
5040



Kelly-Rankin  
Straight  
16cm/6"  
**HHS**  
5041



Kelly-Rankin  
Curved  
16cm/6"  
**HHS**  
5042

**Skin Hooks:** For hooking skin during surgery

**Retractors:** For holding back the lips and cheeks.



25 x 10mm  
22 x 10mm

Farabeuf  
12.5cm  
**4727**

28 x 12mm  
32 x 12mm



23 x 16mm  
20 x 16mm

Farabeuf  
15cm  
**4731**

24 x 16mm  
38 x 16mm




29mm  
28mm  
22mm

Roux  
16.5cm  
**5132**

28mm  
38mm  
45mm

**Retractors:** For holding back the lips and cheeks, and tongue.

**Mouth Gag:** To prop open the mouth and prevent it from closing.





**Retractors:** For holding back the lips and cheeks, and tongue.



**Scalpel Handles:** For holding a mess blade.

**S-42**

|                                                                                      |                                |
|--------------------------------------------------------------------------------------|--------------------------------|
|    | 3<br>12cm<br>5200              |
|    | 3 G<br>12cm<br>5201            |
|    | 4<br>13.5cm<br>5202            |
|    | 5<br>14.5cm<br>5203            |
|  | 5 A<br>14.5cm<br>5204          |
|  | Double Head<br>14.5cm<br>5205  |
|  | 7<br>16.5cm<br>5206            |
|   | 3L<br>21cm<br>5207             |
|   | Scalpel<br>15.2cm<br>5208      |
|   | Scalpel 180°<br>16.2cm<br>5209 |
|   |                                |
|   |                                |
|  |                                |



These tools are used for precision section of surgical sites.

Coupland Handle  
25cm/9.9"  
SASPCH  
5249



Actual Tip Size (Internal diameter)

|              |   |                |       |
|--------------|---|----------------|-------|
| 5250-SASPCT1 | ○ | 1 Coupland Tip | 1.5mm |
| 5251-SASPCT2 | ○ | 2 Coupland Tip | 2.0mm |
| 5252-SASPCT3 | ○ | 3 Coupland Tip | 2.5mm |
| 5253-SASPCT4 | ○ | 4 Coupland Tip | 3.0mm |

Actual Tip Size (Internal diameter)

|                |   |                           |       |
|----------------|---|---------------------------|-------|
| 5254-SASPCT1SH | ○ | 1 Coupland Tip/Side Holes | 1.5mm |
| 5255-SASPCT2SH | ○ | 2 Coupland Tip/Side Holes | 2.0mm |
| 5256-SASPCT3SH | ○ | 3 Coupland Tip/Side Holes | 2.5mm |
| 5257-SASPCT4SH | ○ | 4 Coupland Tip/Side Holes | 3.0mm |

Shown 50% Size.

**BYRD**

Byrd  
Self-Cleaning  
20cm/8"

5258-SASPB

Actual Tip Size (Internal diameter)

|             |   |         |       |
|-------------|---|---------|-------|
| 5259-SASPB1 | ○ | B1 Byrd | 1.5mm |
| 5260-SASPB2 | ○ | B2 Byrd | 2.0mm |
| 5261-SASPB3 | ○ | B3 Byrd | 2.5mm |

Shown 50% Size.

Byrd  
Self-Cleaning  
20cm/8"

5262-SASPCG3


Shown 50% Size.

Actual Tip Size (Internal diameter)

|              |   |            |       |
|--------------|---|------------|-------|
| 5263-SASPF6  | ○ | 6 Frazier  | 1.5mm |
| 5263-SASPF8  | ○ | 8 Frazier  | 2.0mm |
| 5263-SASPF10 | ○ | 10 Frazier | 2.5mm |
| 5263-SASPF12 | ○ | 12 Frazier | 3.0mm |

**Frazier  
Suction Tube**

**Used in Oral Surgery**



Shown 50% Size.

SASP1 SASP2 SASP3 SASP4  
5267A 5267B 5267C 5267D

Oral Surgery  
22cm/8.5"  
1mm

Yandauer  
26.5cm/10.5"

5268-SASPYD

Shown 50% Size.

Andrew-Pychon  
23cm/9"

5269-SASPAP

Shown 50% Size.

Locking Pliers  
16.5cm/6.5"

5270



**Syringes:** An instrument for injecting aesthetics for a cartridge.

**3-Way Syringes:** Supplies air, water, or a mist through separate interval channels created by combining pressurized air with water flow.



Intraligamentary  
Syringe  
1.8ml  
5320




Syringe  
1.8ml  
5322



Syringe  
1.8ml  
5323



Syringe  
1.8ml  
5325



Paroject  
5328



Paroject  
5329



Miscellaneous tools for laboratory.

**S-45**



1  
Lwanson  
Spring Caliper

**CLP1**  
5380



2  
Lwanson  
Spring Caliper

**CLP2**  
5381



Tofflemire  
Junior

**TJ**  
5383



Tofflemire  
Senior

**TS**  
5384



Tofflemire  
Universal

**TU**  
5385



Boley  
Gauge

**BG**  
5386

Our unique tip design makes it mold to the individual shape of your patient's tooth. The smooth handle with opening, allows for more through and easier cleaning. The newly styled handle allows a more uniform transmission while applying pressure, pulling and for rotary movements. More comfortable design can be used by right or left-handed practitioners! Our handles are lighter and are easier to use.



F1  
UPPER  
PREMOLARS

5430



F2  
LOWER  
PREMOLARS

5431



F3  
UPPER  
ANTERIOR

5432



F4N  
UPPER  
ANTERIOR

5433



F6  
UPPER  
3rd MOLAR

5434



Lower Molar  
Modify  
Cow-Horn  
Bead

5430



46XL  
Lower Root  
Extra Long

5436



51XL  
UPPER ROOT

5437



222  
Lower  
Windsoms

5438



301  
Lower  
Roots

5439



300  
Upper Roots

5440

**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



5  
Lower Third  
Molars

5441



15  
Lower  
Molars

5442



10H  
UPPER MOLARS  
HOOK HANDLE

5447



18L  
LOWER MOLARS  
LEFT SIDE

5448



16S  
Lower Molars

5443



17  
Lower Molars

5444



6  
Splitting Forceps

5449



1  
Upper Centrals  
Laterals & Cuspids

5450



16  
Lower Molars  
Cow-Horn Beak

5445



1AF  
Upper  
Anterior

5446

**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



99C  
Upper  
Premolars & Incisors

5451



88C  
Upper Molars  
Left Side

5452



88R  
Upper Molars  
Right Side

5453



10S  
Upper Molars  
Straight-Handle

5454



53L  
Upper Molars  
Left Side

5455



53R  
Upper Molars  
Right Side

5456



23  
Lower Molars  
Cow-Horn Beak

5457



18R  
Lower Molars  
Right Side

5458



69  
Fragments  
& Small Roots  
Upper & Lower

5459



65  
Upper Incisors  
& Roots

5459



Universal  
for Upper  
Teeth & Roots

5461



32  
Upper Molars  
& Premolars  
Strong Beak

5462



**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



151A  
Lower Premolars  
Incisors & Roots

5463



151XAS Children  
Lower Premolars  
Incisors & Roots Canines

5464



151AS  
Split Lower Premolars  
Incisors & Roots Canines

5465



151AF Apical  
Beak Lower Premolars  
Incisors & Roots Canines

5466



101  
Universal  
for Premolars

5467



150  
Upper Premolars  
Incisors & Roots Canines

5468



150AF  
Apical Beak  
Upper Premolars  
Incisors & Roots Canines

5469



150AS  
Split Beak  
Upper Premolars  
Incisors & Roots Canines

5470



150XAS  
Children Split  
Upper Premolars  
Incisors & Roots Canines

5471



150A  
Upper Premolars  
Incisors & Roots Canines

5472



150S  
Upper Premolars  
Incisors & Roots for Children

5473



151  
Upper Premolars  
Incisors & Roots

5474

**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



151S  
Lower Premolars  
Incisors & Roots for Children

5475



210H  
Upper Wisdoms  
With Hook Handle

5476



210S  
Upper Wisdoms  
With Straight Handle

5477



222  
Lower Wisdoms

5478



222AF  
Apical Beak  
Lower Wisdoms

5479



286  
Upper Roots  
Premolars & Incisors

5480

**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



30  
Upper Roots  
& Incisors  
5481



22  
Lower Premolars  
Hawks Bill  
5482



21  
Lower Molars  
5483



20  
Lower Wisdoms  
5484



8  
Lower Premolars  
5485



7  
Lower Premolars  
5486



3  
Upper Incisors  
& Canines  
5487



1  
Upper Incisors  
& Canines  
5488



19  
Upper Wisdoms  
5489



18A  
Upper Molars  
5490



18  
Upper Molars (Left)  
5491



17  
Upper Molars (Right)  
5492



**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.





**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



95  
Upper Molars Left  
5505



121  
Upper Wisdoms  
5506



167 Universal  
Patter for  
Lower Teeth  
Molars & Premolars  
5507



168  
Upper Canines  
5508



74N  
Lower Incisors  
& Roots Canines  
5509



74NF  
Lower Incisors  
& Roots  
5510



79  
Lower Wisdoms  
5511



79A  
Lower Wisdoms  
5512



86C  
Lower Molars  
5513



89  
Upper Molars Right  
5514



90  
Upper Molars Left  
5515



94  
Upper Molars Right  
5516

**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



F - 1  
Upper Universal  
5517



F - 2  
Lower Universal  
5518



F - 3  
Upper Anterior  
5519



F - 4N  
Lower Anterior  
5520



169  
Lower Canines  
5521



233  
Lower Roots  
Parallel Beaks  
5522



235  
Upper Roots  
Parallel Beaks  
5523



MD-1  
Universal Upper Wisdom  
Premolars & Rootis  
5524



MD-4  
Lower Premolars  
5525



22L  
Lower Wisdoms Left  
5526




22R  
Lower Wisdoms Right  
5527

**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



F - 7C  
Lower Molar Modify  
Cow-Horn Beak

5528



F - 8  
Upper Universal

5529



F - 9  
Lower Universal

5530



F - 10  
Lower Molar

5531



F - 11  
Upper Molar Left

5532



F - 12  
Upper Molar Right

5533



F - 13  
Upper Molar Right

5534



F - 14  
Upper Molar Left

5535



F - 4  
Lower Anterior

5536



F - 5  
Lower Anterior

5537



F - 6  
Upper 3rd Molars

5538



F - 7  
Lower 3rd Molars

5539



**Forceps:** For extraction teeth from alveolar bone. American Pattern/English Pattern/Unique Design/ Children's Pattern.

**Root Splinter And Extraction Tweezer:** For removing the residual part of tooth root.

**Sponge & Dressing/Hemostates:** For grasping off blood vessels and removing small piece of root tips.



F - 15  
Upper Premolars  
Universal  
5540



F - 16  
Upper Roots  
5541



F - 17  
Lower Roots  
5542



F - 18  
Upper Roots  
Fragment Narrow  
5543



F - 18TN  
Upper Roots  
Fragment Narrow  
5544



F - 19  
Lower Roots  
Fragment Narrow  
5545



F - 19NW  
Lower Roots  
Fragment Narrow  
5546



46XL  
Lower Roots  
Extra Long  
5547



F - 20  
Upper Premolars  
1ST & 2RD Molars  
5548



F - 49  
Upper Roots  
5549



51XL  
Upper Roots  
5550





F - 1  
Upper Universal  
Pedo

5551



F - 2  
Lower Universal  
Pedo

5552



F - 3  
Upper Anterior  
Pedo

5553



F - 4  
Lower Anterior  
Pedo

5554



F - 5  
Lower Anterior  
Pedo

5555



F - 6  
Upper 3rd  
Pedo

5556



F - 6  
Lower 3rd  
Pedo

5557



F - 16  
Lower Root  
Pedo

5558



F - 17  
Lower Root  
Pedo

5559



B 23  
Lower Molar  
Pedo

5560



150 - SK  
Upper Universal

5561



151 - SK  
Lower Universal

5562

**FRIEDMAN**

**BLUMENTHAL**

**CLEVELAND**



Friedman  
14cm/5-12"  
RF

6000



Blumenthal  
Standard 30°  
14.5cm/5 3/8"  
RBIS

6003



1A  
16.5cm/6 1/2"  
R1A

6006



3H Dean  
17cm/6 3/4"  
R3HD

608



4 Cleveland  
14cm/5 1/2"  
R4C

610



Mini Friedman  
RMF

6001



Blumenthal 45°  
14.5cm/4 3/8"  
RBL45

6004



4A Cleveland  
17cm/6 3/8"  
R4AC

6007



5 Cleveland  
Side Cutting  
17cm/6 3/8"  
R5C

612



5S Cleveland  
Side Cutting  
15cm/6"  
R5SC

613



Mini Friedman 90°  
14.5cm/5 1/2"  
RMF90

6002



Blumenthal 90°  
14.5cm/5 3/8"  
RBL90

6005



Kramer-Nevins  
14.5cm/5 3/8"  
RKN

614



Kerrison Rongeur 3mm  
for widening the sinus window  
RK3mm

6020



Kerrison Rongeur 4mm  
for widening the sinus window  
RK4mm

6022

### DOUBLE ACTION HINGED RONGEURS



18m - (Double Action)  
Side Cutting  
17mm (6-7/8")  
R18MDA

6025



Beyer - (Double Action)  
18cm (7")  
RBEYDA

6030



Boehler - (Double Action)  
15cm (6")  
RBOEDA15

6035



Boehler - (Double Action)  
16cm (6 1/8")  
RBOEDA16

6040



+49 1590 6127425



[info@surgiwelt.com](mailto:info@surgiwelt.com)