

# 2025 VERSION

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments.

[www.surgiwelt.com](http://www.surgiwelt.com)



**Blackline**  
Instruments

  
**SurgiWELT**  
*Medical Instruments Manufacturer*



+49 1590 6127425



[info@surgiwelt.com](mailto:info@surgiwelt.com)



## **COMPANY PROFILE**

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments. With a focus on precision, quality, and innovation.

We provide essential tools for healthcare professionals worldwide. Our products are designed to meet the rigorous demands of modern surgery, ensuring optimal performance and patient safety.

Key areas of expertise:

Manufacturing a wide range of surgical instruments

Commitment to quality and precision

Innovation in surgical instrument design



2  
DEXD2  
BL1134

3CH  
DEXD3CH  
BL1140

5  
DEXD5  
BL1141

TU17/23  
Wilkins/Tufts  
DEXTU17/23  
BL1138

11/12  
Old Dominion University  
DEXD11/12ODU  
BL1147


23/CP-12  
DEP23/CP12  
DEXD11/12ODU  
BL1295



## EXPLORER



## EXCAVATORS





D11TS  
Thinner 21mm  
ERCSD11TS  
BL 1550



Ma57  
Heat Carrier 27mm  
RCSMA57  
BL 1555



GP1-20mm  
RCSGP1  
BL 1544



D11T  
Thinner 22mm  
RCSD11T  
BL 1549



D11  
Standard 23mm  
RCSD11  
BL 1547

Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some case cotton wool rolls. These are different designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

### MOLARS - GENERAL




### UPPER MOLARS



Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some case cotton wool rolls. These are different designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.




Rubber dam clamp applied to grip the cervical area of a tooth to retain rubber dam and in some case cotton wool rolls. These are diffident and designs to fit molars, premolar and incisor.

Martensite stainless steel, Clamps shown 75% actual.

### PRE-MOLARS



### ANTERIOR



Rubber Dam  
Clamp  
Organizing  
Board  
RDCOB  
BL 1877



Rompen Plugger, Mesial/Distal  
BLRPM/D  
BL 1900



Rompen Plugger, Buccal/Lingual  
BLRPB/L  
BL 1901



Jovanovic Plugger  
BLJP  
BL 1902



Labanca Plugger  
BLLBP  
BL 1903



Grisdale Bone Condenser  
BLGBP  
BL 1904



Sinus Curettes 1  
ISC1  
BL 1915



Sinus Curettes 2  
ISC2  
BL 1916



Kramer-Navins  
Obtuse  
IKNO  
BL 1917



Kramer-Navins  
Acute  
IKNA  
BL 1918



Sinus Lift 300  
ISL300  
BL 1919



1/2 Buser Modified  
Periodontal Chisel (4mm/5mm)  
I1/2BMPBC  
BL 1930



2/3 Buser Modified  
Periodontal Chisel (5mm/6mm)  
I2/3BMPBC  
BL 1931



1/3 Buser Modified  
Periodontal Chisel (4mm/5mm)  
I1/3BMPBC  
BL 1932



1/3 Buser Modified  
Periodontal Chisel (4mm/6mm)  
I1/3BMPBC  
BL 1932



3/4 Solt  
ISOLT3/4  
BL 1934



15/16  
Kirkland Knife  
BLKK15/16  
BL 1940



Allen End-Cutting  
Intrасulcular Knife  
BLKAECI  
BL 1941



Allen End-Cutting  
Intrасulcular Knife  
BLKAECI  
BL 1941



1/2 Allen Orban Knife  
Round  
BLKAOR1/2  
BL 1942



1/2 Allen Orban Knife  
BLKAO1/2  
BL 1943



11 Goldman-Fox Knife  
BLG11GFK  
BL 1944



Tunneling Knife 1  
BLTK1  
BL 1950



Tunneling Knife 2  
BLTK2  
BL 1951



Giles De Quincey  
Tunneling Knife 1  
BLGDQTK1  
BL 1952



Giles De Quincey  
Tunneling Knife 2  
BLGDQTK2  
BL 1953



1  
RTNCIGFT1  
BL 3230



6 (Flexible Thin)  
RTNCIGFT6  
BL 3235



2  
RTNCIGFT2  
BL 3231



Micro-Mini  
RTNCIPCS  
BL 3236



3 (Extra Flexible)  
RTNCIGFT3  
BL 3232



Mini 1  
RTNCIGFTM1  
BL 3237



4 (Extra-Flexible)  
RTNCIGFT4  
BL 3233



Mini 3 (Extra-Thin)  
RTNCIGFTM3  
BL 3238



5 (Flexible-Thin)  
RTNCIGFT5  
BL 3234



Mini 4 (Extra-Thin)  
RTNCIGFTM4  
BL 3239



Ab1  
Boghosian  
RTNPFIAB1  
BL 3240



Ab2  
Boghosian  
RTNPFIAB2  
BL 3241



Interproximal Carver  
RTNCVIPC  
BL 3242



W3  
RTNPFIW3  
BL 3243



A6  
RTNPPIA6  
BL 3244



4/5 Gregg  
RTNPFIG4/5  
BL 3245



Contouring Instruments  
RTNCFIR/L  
BL 3246



8A  
RTNPFI8A  
BL 3247

### BURNISHERS



21B  
RTNBB21B  
BL 3250



27/29  
RTNBB2729  
BL 3251



2 Ladmore  
RTNBBL2  
BL 3252



3 Ladmore  
RTNBBL3  
BL 3253



Small/Medium  
Ball Burnisher  
RTNBBS/M  
BL 3254

### FREEMAN BURNISHERS



Freedman (Duckhead)  
RTNPCCI  
BL 3255



Freedman (Small)  
Contact Forming  
RTNFCIS  
BL 3256



Freedman (Large)  
Contact Forming  
RTNFCIL  
BL 3257



A  
Cosmetic Contouring  
RTNCCIA  
BL 3260



B  
Cosmetic Contouring  
RTNCCIB  
BL 3261



F  
Cosmetic Contouring  
RTNCCIF  
BL 3262



C  
Cosmetic Contouring  
RTNCCIC  
BL 3263



G  
Ridge & Embrasure  
Shaping Instruments  
RTNCCIG  
BL 3264



D  
Cosmetic Contouring  
RTNCCID  
BL 3265



H  
Occlusal Anatomy Instruments  
RTNCCIH  
BL 3266



E  
Cosmetic Contouring  
RTNCCIE  
BL 3267



I  
Cosmetic  
Packing Instruments  
RTNCCII  
BL 3268

### ANTERIOR KIT



3  
EXTRA-FLEX  
RTNCIGFT3  
BL 3270



Mini 1  
RTNCIGFTM1  
BL 3272



Micro-Mini  
RTNCIPCS  
BL 3274



Medium Placing/Condensing  
RTNCIPCM  
BL 3276



Large  
Placing/Condensing  
RTNCIPCL  
BL 3278

### POSTERIOR KIT



OT  
Tanner  
RTNPLGOT  
BL 3271



3  
Hollenback  
RTNPLGH3  
BL 3273



5A  
RTNPLG5A  
BL 3275



Small/Medium  
Contact Forming  
RTNCFIS/M  
BL 3277



Medium/Large  
Contact Forming  
RTNCFIM/L  
BL 3279



Probe 0.2mm  
BLPT23  
BL 3582



Heidemann Spatula  
BLHST2  
BL 3580



Probe 0.2mm  
BLPT11  
BL 3581



Plugger (Willand Shape)  
Q.4mm  
BLPWST11  
BL 3585



Spatula T1  
BLST1  
BL 3586



Spatula (Cobra Shape)  
0.4mm  
BLSCST2  
BL 3587



Spatula T4  
BLST4  
BL 3588



Plugger (Straight)  
0.5-0.8mm  
RPS  
BL 3600



Plugger T1  
Bent Shape  
0.5-0.8mm  
RPBST1  
BL 3601



Plugger T2  
Bent Shape  
1.0-1.3mm  
RPBST2  
BL 3602



Plugger T3  
Bent Shape  
0.5-0.8mm  
RPBST3  
BL 3603



Plugger T4  
Bent Shape  
1.0-1.3mm  
RPBST4  
BL 3604



Plugger T5  
Columbia Shape  
0.6mm  
RPCST5  
BL 3610



Plugger T6  
Columbia Shape  
0.8mm  
RPCST6  
BL 3611



Plugger T12  
Wiland Shape  
0.6mm  
RPWST12  
BL 3612



Plugger T13  
Wiland Shape  
0.8mm  
RPWST13  
BL 3613



Plugger T14  
Wiland Shape  
1.0mm  
RPWST14  
BL 3614



Excavator T1  
Wiland Shape  
0.8mm  
REWST1  
BL 3620



Excavator T2  
Wiland Shape  
1.0mm  
REWST2  
BL 3621



Excavator T3  
Wiland Shape  
1.3mm  
REWST3  
BL 3622



Excavator T4  
Columbia Shape  
0.8mm  
RECST4  
BL 3623



Excavator T5  
Columbia Shape  
1.0mm  
RECST5  
BL 3624



Excavator T6  
Columbia Shape-1.3mm  
RECST6  
BL 3625



Excavator T61/2  
Straight/1.0mm  
REST61/2  
BL 3626



Excavator T63/4  
Straight/1.3mm  
REST63/4  
BL 3627



Curette 13/14  
Columbia Shape  
RCCS13/14  
BL 3628



Scaler 13/14  
Columbia Shape  
RCCS13/14  
BL 3629



Buser Periosteal Elevator  
BLBPE  
BL 3650



Allen Periosteal Elevator  
BLAPE  
BL 3651



Allen Periosteal Elevator, Anterior  
BLAPEA  
BL 3652



9 Molt Periosteal Elevator  
BL9MPE  
BL 3653



3 Prichard Periosteal Elevator  
BL3PPE  
BL 3654



2/4 Molt Curette  
BLMC2/4  
BL 3655

## PERIOSTEAL



Swiss Perio  
Modified PS24GSP Periosteal  
BLSPMP24GSPP  
BL 3670



Swiss Perio  
Modified P24GSP Periosteal  
BLSPMP24GSICP  
BL 3671



Swiss Periosteal  
Modified P24GSP Periosteal2  
BLSPMP24GSICPP2  
BL 3672



Swiss Periosteal  
Modified P24GSP Periosteal2  
BLSPMP24GSICPP2  
BL 3672



Swiss Periosteal  
Modified Schlee Periosteal  
BLSPMSP  
BL 3673



Modified Labanca Periosteal  
BLMPL  
BL 3675



Urban Mini Me Periosteal  
PLUMMP  
BL 3680



Buser/Modified PH2MBH (Large)  
BLBMPH2MBHL  
BL 3682



Buser/Modified PH2MBH (Large)  
BLBMPH2MBHL  
BL 3682




152-Kramer-Nevius Periosteal  
BLKN152P  
BL 3683



24G Periosteal  
BL24GP  
BL 3684



14 Goldsman-Fox Periosteal  
BL14GFP  
BL 3685



**EZ-LUXATING ELEVATORS**

**SYNDESMOTOMES**



EX-Luxating  
Straight, Spade  
BLEZLSS  
BL 3700



EZ-Luxating  
Small, Straight  
BLEZLSS  
BL 3701



EZ-Luxating  
Small, Distal  
BLEZLSD  
BL 3703



EZ-Luxating  
Large, Distal  
BLEZLSD  
BL 3705



6 Right Bernard  
Syndesmotosomes  
BL6RBS  
BL 3710



EZ-Luxating  
Small, Curved  
BLEZLSC  
BL 3702



EZ-Luxating  
Small, Mesial  
BLEZLSM  
BL 3704



EZ-Luxating  
Large, Mesial  
BLEZLLM  
BL 3706



7 Left Bernard  
Syndesmotosomes  
BL7LBS  
BL 3711



11 Straight Barnard  
Syndesmotosomes  
BL11SBS  
BL 3712

Weider Tongue  
Small

**BL 5150**

**Weider Tongue  
Retractor, Small  
BLWTRS  
BL 5150**

Weider Tongue  
Large

**BL 5151**

**Weider Tongue  
Retractor, Large  
BLWTRL  
BL 5151**

Bishop

**BL 5140**

**BLBIS  
BL 5140**

University of  
Minnesota

**BL 5133**

**BLUOM  
BL5133**

Lindo  
Serrated

**BL 5134**

**BLLS  
BL5134**

Front Round

**BL 5135**

**BLFR  
BL5135**

Front Cut

**BL 5136**

**BLFC  
BL5136**

Cawood  
Minnesota

**BL 5137**

**BLCWM  
BL5137**





Luxator, Curved  
3mm  
BLLXC3  
BL 4006



Luxator, Straight  
4mm  
BLLXS4  
BL 4007



Luxator, Curved  
4mm  
BLLXC4  
BL 4008



Luxator, Straight  
5mm  
BLLXS5  
BL 4009



Luxator, Curved  
5mm  
BLLXC5  
BL 4010



Luxator, Straight  
2mm  
BLLXS2  
BL 4000



Luxator, Curved  
2mm  
BLLXC2  
BL 4002



Luxator, Straight  
2.5mm  
BLLXS2.5  
BL 4003



Luxator, Curved  
2.5mm  
BLLXC2.5  
BL 4004



Luxator, Straight  
3mm  
BLLXS3  
BL 4005

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction

**B-27**

## HEIDBRINK

## APICAL



1  
Heidbrink  
**EHB1**  
BL 4316



2  
Heidbrink  
**EHB2**  
BL 4317



3  
Heidbrink  
**EHB3**  
BL 4318



13/14  
Heidbrink  
**EHB13/14**  
BL 4312



9R  
**E9R**  
BL 4323



9L  
**E9L**  
BL 4324

These instruments are used to loosen the tooth from the periodontal ligament and ease extraction

**B-28**

## DAVIS



10  
Davis  
Double-End  
**ED10**  
BL 4343

11  
Davis  
Double-End  
**ED11**  
BL 4344

## HOWARD



2  
Howard  
**EHO2**  
BL 4332

3  
Howard  
**EHO3**  
BL 4333





Periotome,  
Posterior (DE)  
BLPPDE  
BL 4347

Periotome,  
Posterior, Right (SE)  
BLPPRSE  
BL 4348

Periotome,  
Posterior, Left (SE)  
BLPLPSE  
BL 4349

Periotome,  
Anterior (DE)  
BLPADE  
BL 4350

Periotome,  
Anterior, Angled (SE)  
BLPAASE  
BL 4351

Periotome,  
Anterior, Straight (SE)  
BLPASDE  
BL 4352



 +49 1590 6127425

 [info@surgiwelt.com](mailto:info@surgiwelt.com)