

2025 VERSION

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments.

www.surgiwelt.com



Periodontics
Instruments


SurgiWELT
Medical Instruments Manufacturer



+49 1590 6127425



info@surgiwelt.com



COMPANY PROFILE

Surgiwelt is a dedicated manufacturer specializing in the creation of surgical instruments. With a focus on precision, quality, and innovation.

We provide essential tools for healthcare professionals worldwide. Our products are designed to meet the rigorous demands of modern surgery, ensuring optimal performance and patient safety.

Key areas of expertise:

Manufacturing a wide range of surgical instruments

Commitment to quality and precision

Innovation in surgical instrument design



INDEX

DOUBLE ENDED



9.5mm
(3/8")

9.5mm
(3/8")

9.5mm
(3/8")

SINGLE ENDED



9.5mm
(3/8")

9.5mm
(3/8")

9.5mm
(3/8")

Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.



H6/H7
Hygienist
SH6/7
2630



6/7
Goldman
SHG6/7
2631



H5/U15
Hygienist-Towner
SH5/U15
2632



U15/30
Towner-Jacquette
SU15/30
2633



H5/33
Hygienist-Jacquette
SH5/33
2634



U15/33
Towner-Jacquette
SU15/33
2635



U15
Towner
SU15
2636

Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.



Nevi1
Anterior

SCNEVI19
2637



Nevi 1/H5
Anterior

SCNEV1/H56
2638



Nevi 2
Posterior

SCNEVI29
2639



Nevi 3
Posterior

SCNEVI39
2640



Nevi 4
Posterior

SCNEVI49
2641

Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.



30/33
Jacquette
SJ30/33
2642



31/32
Jacquette
SJ31/32
2643



34/35
Jacquette
SJ34/35
2644



H3/H4
Jacquette
SH3/4
2645



14/15
Posterior
Jacquette
SJ14/15
2646

Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.



204

S204
2647



204S

S204S
2648



204SD

S204SD
2649



204IUFW

SIUFW204
2650



1
Springfield
Anterior

SSP1
2651



2/3
Springfield

SSP2/3
2652

Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.

P-7

PROSTERIOR

ANTERIOR



Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.

ANTERIOR



Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.

P-12

SICKLE SCALERS

BALL SCALERS

SRANE-KAPLAN



152
Cumine

PSSCM152
2694



CI-2/3

PSSCI-2/3
2695



CI-1

PSSCI-1
2696



B-1

PSSB-1
2697



B-2

PSSB
2698



B-5/6

PSSB-5/6
2699



B-2/3

PSSB-2/3
2700



CK-6

PSSCK-6
2701



CK-6D

PSSCK-6D
2702



Sickle Scalers: Sickle Scalers are strong instruments that remove large calculus deposits. A sickle scaler 2 cutting edges per working end. Sickle scalers may be used on all aspects of the tooth crown, facial, lingual, mesial and distal. The working end is triangular in cross-section. Each of the lateral surfaces meets the face at an internal angle of 70° to 80°.

P-13



11/12
Scaler

PSS11/12
2705

11A/12A
Scaler

PSS11/12A
2706

11/12S
Modified
McCall

PSS11/12S
2707

13S/14S
Pointed
Curette

PSC13/14S
2708

MTC 13/14
Pointed
Curette

PSCTC13/14
2709

17S/18S
Curette

PSC17/18
2710

19/20
Curette

PSC19/20
2711

One side is sickle scalar and the other is a universal curette

P-14

Goldman-Fox Set



One side is sickle scalar and the other is a universal curette

P-15



GF-4
Curette

PSGF4C

2766



GF-21S/33
Comination

PSGF21S/33C

2767


Cattoni Set



103/106

PS103/106

2768



107/108

PS107/108

2769



104/105

PS104/105HC

2770



107/108S

PS107/108SM

2771

Blades are sharpened on both sides. Blade curved at 90° to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

P-16

The langer cures combine the shank design of the gracey cures with the 90° cutting edge of the universal cures, edged on both sides. As the langer curette 1-2 is applicable for the treatment of lower bicuspid and molars, and langer curette 3-4 for upper bicuspid and molars, as well as the langer curette 5-6 for upper and lower front teeth, the model 17-18 is specially adapted for all surfaces in the area of upper and lower jaw, specially at small stoma.



1/2
Langer
(Shown)
PUCL1/2
2820

1/2
Langer
After 5
PUCL1/2A5
2822

1/2
Langer
Rigid
PUCL1/2R
2821

1/2
Langer
Mini 5
PUCSL1/2M5
2823



3/4
Langer
(Shown)
PUCL3/4
2824

3/4
Langer
After 5
PUCL3/4A5
2826

3/4
Langer
Rigid
PUCL3/4R
2825

3/4
Langer
Mini 5
PUCSL3/4M5
2827



5/6
Langer
(Shown)
PUCL5/6
2828

5/6
Langer
Rigid
PUCL5/6R
2829

5/6
Langer
After 5
PUCL5/6A5
2830

5/6
Langer
Mini 5
PUCSL5/6M5
2831



17/18
Langer
(Shown)
PUCL17/18
2832

17/18
Langer
Rigid
PUCL17/18R
2833

Blades are sharpened on both sides. Blade curved at 90° to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

P-18

Anterior



IU13/14
Indiana
University
PUCIU13/14
2834

Posterior



IU17/18
Indiana
University
PUCIU17/18
2835

Curette-Chisel Combination



1/2
Bunting
PCCCB1/2
2836



5/6
Bunting
PUCBN5/6
2837



YG7/8
Younger Good
PUCYG7-8
2838

Blades are sharpened on both sides. Blade curved at 90° to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

P-19

Barnhart



1/2
Barnhart
(Shown)
PUCBH1/2
2839

1/2
Barnhart
Rigid
PUCBH1/2R
2840



5/6
Barnhart
(Shown)
PUCBH5/6
2841

5/6
Langer
Mini 5
PUCBH5/6R
2842

Columbia



2R/2L
Columbia
PUCC2R/2L
2843



4R/4L
Columbia
PUCC4R/4L
2844



13/14
Columbia
Columbia
PUCC13/14
2845

13/14
Columbia
Rigid
PUCC13/14R
2846

13/14
Columbia
Curved Tip
PUCC13/14C
2847

Blades are sharpened on both sides. Blade curved at 90° to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

P-20



13S/14S

Bates

PUCBT13/14S

2848



HU/12

Hutchinson

PUCHU1/2

2849



3/4

Rule (23-24)

PUCR3/4

2850



UC3/4

Ratcliff

PUCUR3/4

2851



411/412

PUC411/412

2852

Blades are sharpened on both sides. Blade curved at 90° to shank with a rounded toe. Designed so that the working ends can be adopted to all tooth surfaces of all regions of the mouth with one double ended instrument.

P-21



10/11
Loma Linda
(Shown)

PUCLL10/11
2853

10/11
Loma Linda
Rigid

PUCLL10/11R
2854

10/11
Loma Linda
Slim

PUCLL10/11S
2855



413/414

PUC413/414
2856



1/2
Mallery

PUCMAL1/2
2857



31/32
Goldman

PUCGAL31/32
2858



ME-1
[10-11]

PUCME1
2859

Curettes are used supragingivally and subgingivally on enamel and root surfaces. These are two major sub-types of curettes, Universal curettes and area-specific curettes. Universal curettes can be used to remove medium-sized and small calculus deposits. Area-Specific curettes are designed to remove light calculus and endotoxins from root surfaces. Gracey curettes is an area-specific curettes. It has offset blade that provides a perfect working angulation for the teeth surface.



	1/2 (Shown)	3/4 (Shown)	5/6 (Shown)	7/8 (Shown)	9/10 (Shown)		
	PCG1/2	PCG3/4	PCG5/6	PCG7/8	PCG9/10		
	2910	2912	2914	2917	2920		
1/2 Rigid	3/4 Rigid	5/6 Rigid	7/8 Rigid	9/10 Rigid	5/6 Extra Rigid	7/8 Extra Rigid	
PCG1/2R	PCG3/4R	PCG5/6R	PCG7/8R	PCG9/10R	PCG5/6ER	PCG7/8ER	
2911	2913	2915	2918	2921	2916	2919	

Curettes are used supragingivally and subgingivally on enamel and root surfaces. These are two major sub-types of curettes, Universal curettes and area-specific curettes. Universal curettes can be used to remove medium-sized and small calculus deposits. Area-Specific curettes are designed to remove light calculus and endotoxins from root surfaces. Gracey curettes is an area-specific curettes. It has offset blade that provides a perfect working angulation for the teeth surface.



11/12
(Shown)

PCG11/12
2922

13/14
(Shown)

PCG13/14
2925

15/16
(Shown)

PCG15/16
2928

17/18
(Shown)

PCG17/18
2930

11/12
Rigid

PCG11/12R
2923

13/14
Rigid

PCG13/14R
2926

15/16
Rigid

PCG15/16R
2929

17/18
Rigid

PCG17/18R
2931

11/12
Extra Rigid

PCG11/12ER
2924

13/14
Extra Rigid

PCG13/14ER
2927

Rigid Gracey: The rigid gracey has a stronger and less flexible shank than the standard one. So it can be used for heavy calculus removal without using additional scalars.

P-24

Gracey Mesial Curettes

Gracey Prophy Curettes



11/14
Mesial Distal
(Shown)

PGMDC11/14
2932



12/13
Mesial Distal
(Shown)

PGMDC12/13
2934

11/14
Mesial Distal
Rigid

PGMDC11/14R
2933

12/13
Mesial Distal
Rigid

PGMDC12/13R
2935



P 1/2

PGPC1/2
2936



P 3/4

PGPC3/4
2937



P 5/6

PGPC5/6
2938




P 7/8

PGPC7/8
2939



After-5 Curettes: 1. Terminal shank elongated by 3mm for access into deep periodontal pockets and root surfaces of 5mm or more and better clearance around the crown.
2. Blade thinned by 10% to ease gingival insertion and reduce tissue distention.
3. Area-Specific designed patterned after the standard graceys offers exact adaption and control.
4. Available in Rigid and Extra Rigid designs for greater strength when needed.

P-25



1/2

3/4

5/6

7/8

11/12

13/14

15/16

11/14
Mesial Distal

12/13
Mesial Distal

A5G1/2
2940

A5G3/4
2942

A5G5/6
2944

A5G7/8
2947

A5G11/12
2950

A5G13/14
2953

A5G15/16
2956

A5G11/14
2958

A5G12/13
2960

Rigid

Rigid

Rigid

Rigid

Rigid

Rigid

Rigid

Rigid

Rigid

A5G1/2R
2941

A5G3/4R
2943

A5G5/6R
2945

A5G7/8R
2948

A5G11/12R
2951

A5G13/14R
2954

A5G15/16R
2957

A5G11/14R
2959

A5G12/13R
2961

Extra Rigid

Extra Rigid

Extra Rigid

Extra Rigid

A5G5/6ER
2946

A5G7/8ER
2949


A5G11/12ER
2952

A5G13/14ER
2955



After-5 Curettes: 1. Terminal shank elongated by 3mm for access into deep periodontal pockets and root surfaces of 5mm or more and better clearance around the crown.
2. Blade thinned by 10% to ease gingival insertion and reduce tissue distention.
3. Area-Specific designed patterned after the standard graceys offers exact adaption and control.
4. Available in Rigid and Extra Rigid designs for greater strength when needed.

P-26



1/2

M5G1/2
2962

Rigid

M5G1/2R
2963

3/4

M5G3/4
2964

Rigid

M5G3/4R
2965

5/6

M5G5/6
2966

Rigid

M5G5/6R
2967

7/8

M5G7/8
2968

Rigid

M5G7/8R
2969

11/12

M5G11/12
2970

Rigid

M5G11/12R
2971

13/14

M5G13/14
2972

Rigid

M5G13/14R
2973

15/16

M5G15/16
2974

Rigid

M5G15/16R
2975

11/14
Mesial Distal
M5G11/14
2976

12/13
Mesial Distal
M5G12/13
2977



Hoe Scalars: Use to dislodge heavy sepramarginal calculus deposits use a pull stroke.

File Scalars: Use to crush and remove heavy calculus deposits. Use a pull stroke.

P-27

File Scalars



3/7
Orban
Buccal/Lingual
PFS-OBL3/7
3030

6/7
Orban
Buccal/Lingual
PFS-OBL6/7
3031

12/13
Orban
Buccal/Distal
PFS-OMD12/13
3032

5/11
Hirschfeld
Mesial Distal
PFS-HMD5/11
3033

9/10
Hirschfeld
Anterior
PFS-HA9/10
3034

Hoe Scalars: Use to dislodge heavy sepramarginal calculus deposits use a pull stroke.

File Scalars: Use to crush and remove heavy calculus deposits. Use a pull stroke.

P-28

Hoe Scalars



XS-15
Towner

PHSXS15
3035



7K

PHS7K
3036



6/7
Orban
Mesial Distal

PHSO6/7
3037



8/9
Orban
Buccal/Lingual

PHSO8/9
3038



U15/13R

PHSU15/13R
3039

For Incision that remove or raconteur soft tissue. Also useful to excise inter-proximal tissue.

P-29



15K/16K
Kirkland

PK15K/16K
3080



7
Gold-Fox

PKG7
3081



8
Gold-Fox

PKG8
3082



9
Gold-Fox

PKG9
3083



11
Gold-Fox

PKG11
3084



3/4
Buck

PKB3/4
3085



5/6
Buck

PKB5/6
3086

For Incision that remove or raconteur soft tissue. Also useful to excise inter-proximal tissue.

P-30



1/2
Soft

PKSOLT1/2
3087



1/2
Sanders

PKS1/2
3088



1/2
Orban

PKQ1/2
3089



USC
Tow near
19/20

PKUSCT19/20
3090



1/2
Orban
Knife

PKOK1/2
3091



CK-3
Crane
Kaplan

PKCK3
3092



Dr Allen
End-Cutting
Instrument Knife

PKDA
3093





3 Kramer-Nevins
7.5mm
PKN3
3100

55 Kramer-Nevins
5.5mm
PKN55
3101

75 Kramer-Nevins
7.5mm
PKN75
3102

4 Buser
Round Tip-4mm
PBS4
3103

6 Buser
Pointed Tip-6mm
PBS6
3104

1 Ochesenbein
PCO1
3105

2 Ochesenbein
PCO2
3106





+49 1590 6127425



info@surgiwelt.com